MicronJet™

Microneedles for Intradermal Injection

EN	4
DE	8
ES	12
FR	16
IT	20
PT	24
PL	28

NL	32
CS	36
EL	40
SV	44
TR	48
AR	52


Device preparation for injection

- Draw injection solution into the syringe without purging.
- Open the soft blister and **firmly connect** the MicronJet™ to the syringe.
- Remove the plastic cover.
- Purge air and excess fluid leaving the required volume in the syringe.
- Make sure injection site is dry prior to injection.

The MicronJet™ Microneedles for intradermal injection

- **Hold syringe flanges** between thumb and fingers.
- The **blue line** must face you at all times.
- Syringe flanges are not the orientation guide as they may not always be parallel to the blue line; use the blue line for orientation.

- **Stretch the skin** downwards throughout the entire injection.
- Insert the device at a **45°-70° angle** (flatten as required to reduce injection pressure)
- Press the device into the skin until a **shallow indentation** is formed.

- **Inject slowly**, expect strong resistance.
- A **white wheal** (bleb) will form at the injection site.
- Upon completion, **wait 3 seconds** and remove syringe.

MicronJet™ Instructions for Use

Description:

The MicronJet™ is a micro-needle device that is mounted on a standard syringe in the same manner as a conventional needle.

The MicronJet™ comprises an integral array of micro-needles:

- The MicronJet™ 600 – 0.6mm length, 3 micro-pyramids
- MicronJet™ 800³ – 0.8mm length, 3 micro-pyramids.

The MicronJet™ is used to inject liquid substances, allowing for controlled intradermal delivery in any procedure which requires administration of substances to the dermal compartment.

The MicronJet™ is latex free.

Indications:

The MicronJet™ is intended for intradermal injections of any substance or drug approved for delivery by this delivery route.

Contraindications and limitations:

- The MicronJet™ should not be used on skin abrasions, open wounds, cuts and scars.

- The MicronJet™ should not be used on rashes, skin infections or any other area of damaged or diseased skin.

Warnings and precautions:

- The MicronJet™ is intended for single use, by healthcare professionals only.
- The MicronJet™ can be used in all ages and populations.
- Re-use may lead to infection or other illness/injury.
- Used or un-used devices should not be re-capped or removed from syringe unless there is no alternative or such action is required by a specific medical procedure.
- The MicronJet™ is not intended for aspiration of liquids.
- If the device appears to be damaged or broken, discard and use a new device.
- Any serious incidence should be reported to the manufacturer and the competent authority.

Adverse events:

Local adverse reactions, including edema, erythema and discoloration of skin at the injection site have been observed, most likely related to the injected substance, and not to the MicronJet™.

MicronJet™
Microneedles for Intradermal Injection

Adverse reactions related to use of the MicronJet™ include small risk of micro bleeding [self-limiting blood dotting, ceasing within normal clotting time]. Other potential device related adverse events may rarely include local infection at the injection site.

Storage and handling:

- MicronJet™ should be stored at -10 to 30°C [14 to 86 °F].
- The MicronJet™ is individually packed in a soft blister, sterilized by ethylene oxide (EtO).
- Do not use the MicronJet™, in the event that the soft blister is either open and/or damaged.
- Dispose of all devices in accordance with national/institutional guidelines, and in their absence, discard in a sharps container.

Supply method:

The device is supplied as a single sterile MicronJet™ device, packed in a soft blister for single use.

USA only - Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed health care professional.

Symbol definitions

	MANUFACTURER		DO NOT USE IF PACKAGE IS DAMAGED AND CONSULT INSTRUCTIONS FOR USE		MEDICAL DEVICE
	AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY/ EUROPEAN UNION		SINGLE STERILE BARRIER SYSTEM, STERILIZED USING EtO		UNIQUE DEVICE IDENTIFIER
	USE-BY/EXPIRATION DATE		SINGLE STERILE BARRIER SYSTEM WITH PROTECTIVE PACKAGING OUTSIDE, STERILIZED USING EtO		PACKAGING UNIT
	LOT/BATCH NUMBER		KEEP DRY		TEMPERATURE LIMIT
	CATALOGUE NUMBER		DO NOT REUSE		PRESCRIPTION ONLY (USA)
	MODEL TYPE		CONSULT INSTRUCTIONS FOR USE		CE CONFORMITY MARKING
	DATE OF MANUFACTURE		NON PYROGENIC		FRAGILE, HANDLE WITH CARE
	DO NOT RESTERILIZE				KEEP AWAY FROM SUNLIGHT

Manufacturer:
NanoPass Technologies Ltd.
3 Golda Meir St., Nes Ziona
7403648,
Israel
Tel +972-8-9462905
E-mail: info@nanopass.com

EU Representative:
MedNet EC-REP GmbH
Borkstrasse 10-48163, Muenster,
Germany
Phone +49-251-322660
Fax +49-251-3226622

*For any further assistance/support/questions,
please contact the manufacturer.*

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MedNet EC-REP GmbH Borkstrasse 10-48163,